

SIEMENS MAGNETOM Verio syngo MR B17

\\USER\PHYSICS\Pete\Yuqi_qDESS_240225\qDESS_PJL_original_lores

TA: 3:23 PAT: 2 Voxel size: 0.7x0.7x1.5 mm Rel. SNR: 1.00 USER: qDESS_PJL

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slab group 1	
Slabs	1
Dist. factor	20 %
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	80
FoV read	150 mm
FoV phase	100.0 %
Slice thickness	1.50 mm
TR	25.50 ms
TE	6.00 ms
Averages	1
Filter	None
Coil elements	KN

Contrast

Flip angle	20 deg
Fat suppr.	Water excit. fast
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

Resolution

Base resolution	224
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	Off
Slice partial Fourier	Off
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	24
Accel. factor 3D	1
Matrix Coil Mode	Auto (Triple)
Reference scan mode	Integrated
Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off
Raw filter	Off

Elliptical filter

Off

Geometry

Multi-slice mode	Sequential
Series	Interleaved
Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

System

Body	Off
KN	On
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Sagittal
Rotation	0.00 deg
F >> H	150 mm
A >> P	150 mm
R >> L	120 mm

Inline

Subtract	Off
Std-Dev-Sag	Off
Std-Dev-Cor	Off
Std-Dev-Tra	Off
Std-Dev-Time	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Sequence

Introduction	Off
Dimension	3D
Elliptical scanning	On
Asymmetric echo	Off
Bandwidth	199 Hz/Px
Flow comp.	No
Segments	1
RF pulse type	Normal
Gradient mode	Whisper
Excitation	Slab-sel.
qDESS spoiler duration	3.62 ms
qDESS spoiler moment	31.30 mT/m.ms
Balanced readout	On

SIEMENS MAGNETOM Verio syngo MR B17

\\USER\PHYSICS\Pete\Yuqi_qDESS_240225\qDESS_PJL_proposed_lores

TA: 4:00 PAT: 2 Voxel size: 0.7x0.7x1.5 mm Rel. SNR: 1.00 USER: qDESS_PJL

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slab group 1	
Slabs	1
Dist. factor	20 %
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	80
FoV read	150 mm
FoV phase	100.0 %
Slice thickness	1.50 mm
TR	30.14 ms
TE	4.03 ms
Averages	1
Filter	None
Coil elements	KN

Contrast

Flip angle	41 deg
Fat suppr.	Water excit. fast
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

Resolution

Base resolution	224
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	Off
Slice partial Fourier	Off
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	24
Accel. factor 3D	1
Matrix Coil Mode	Auto (Triple)
Reference scan mode	Integrated
Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off
Raw filter	Off

Elliptical filter

Off

Geometry

Multi-slice mode	Sequential
Series	Interleaved
Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

System

Body	Off
KN	On
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Sagittal
Rotation	0.00 deg
F >> H	150 mm
A >> P	150 mm
R >> L	120 mm

Inline

Subtract	Off
Std-Dev-Sag	Off
Std-Dev-Cor	Off
Std-Dev-Tra	Off
Std-Dev-Time	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Sequence

Introduction	Off
Dimension	3D
Elliptical scanning	On
Asymmetric echo	Off
Bandwidth	199 Hz/Px
Flow comp.	No
Segments	1
RF pulse type	Normal
Gradient mode	Fast
Excitation	Slab-sel.
qDESS spoiler duration	3.62 ms
qDESS spoiler moment	31.30 mT/m.ms
Balanced readout	On